



# Safety. Innovation. Freedom.

## PPE DEFENDS Rapid (15 Minutes) Covid-19 Antigen Test

### FEATURES

- Lateral Flow Test
- Time to result = 15 minutes
- Storage Temperature Range 2 - 30° C
- Shelf Life 24 months
- Specification Type nasopharyngeal swab
- Highly accurate 98.74%
- UK government Approved

### REASSURANCE YOU CAN TRUST



### CERTIFICATION



IVD CE notification  
p109656 Edmazzde  
2004 8 0 2020-09-11-  
medical device safety service



International  
Organization for  
Standardization  
(ISO 13485:2016)



MHRA  
CTK Biotech branded COVID-19  
IgG/IgM Rapid Test Device  
product can be placed on the  
UK market for professional use



World Health  
Organization  
CTK Biotech are currently  
recognised supplier by the WHO

## CONTINUITY OF BUSINESSES IN A COVID SAFE ENVIRONMENT



UK IVD notification  
p109656 Edmazzde  
2004 8 0 2020-09-11-  
medical device safety service

International  
Organization for  
Standardization  
(ISO 13485:2016)

MHRA  
CTK Biotech branded COVID-19  
IgG/IgM Rapid Test Device  
product can be placed on the  
UK market for professional use

World Health  
Organization  
CTK Biotech are currently  
recognised supplier by the WHO

Transport Companies, Offices, Factories, Schools, Food Processing, Hospitality and many other Businesses can deploy **Rapid Testing** as part of their COVID-19 safety and security procedures.

Providing another layer of reassurance, It will enhance confidence in the business from employees, customers and suppliers alike.

**PPE Defends** supplies the Healgen COVID-19 Antigen Test (Swab). It is intended to aid in the rapid diagnosis of SARS-CoV-2 infections.

Unlike many rapid test products, the Healgen test has been put through testing at Porton Down as part of the UK Governments assessment program, and as a result is being used in the Liverpool mass testing project, which has been performing more than 12,000 tests per day.

[www.ppedefends.com](http://www.ppedefends.com)

Tel: 01276 984225

[sales@ppedefends.com](mailto:sales@ppedefends.com)